



固位钉(膜钉TACK) 固位螺丝(钛钉SCREW)




图片仅供参考,产品以实物为准
最终解释权归商家所有

产品特性与优势



产品特点:



适应症

为骨增量提供辅助支撑、固定作用



作用

可以有效固定屏障膜/钛膜



材质

选用优质钛合金材质生物相容性良好



特性

规格丰富,可辅助钛膜多点位固定



便捷

简化手术操作
统一灭菌交付,即开即用

尺寸规格表


固位螺丝 (钛钉SCREW)



固位螺丝 (Screw)	型号 Model	钉直径(mm) D	长度(mm) L	帽直径(mm) Ø
	LDGL1404	1.4	4	2.66
	LDGL1406	1.4	6	2.66
	LDGL1408	1.4	8	2.66
	LDGL1410	1.4	10	2.66
	LDGL1604	1.6	4	2.66
	LDGL1606	1.6	6	2.66
	LDGL1608	1.6	8	2.66
	LDGL1610	1.6	10	2.66

固位钉 (膜钉TACK)



固位钉 (Tack)	型号 Model	钉直径(mm) D	长度(mm) L	帽直径(mm) Ø
	LDMD1003	0.7	3.0	2.5
	LDMD1004	0.7	4.0	2.5

工具型号/规格

固位钉持钉器



70mm



92mm



66.5mm

157mm



74mm

164mm

型号

BTA01

BTA02

BTA03

BTA04

固位螺丝/螺丝刀



20mm



26mm



60mm

型号

BSA01

BSA02

BSA03

工具型号/规格

定位钻



型号

GPD08 GPD12 GPD14

滋养孔钻



型号

N104 N124 N144

牙骨锤



型号

800/510

蓝色手柄



型号

BHD01

屿枫® 膜钉TACK



套餐A



一次购买5盒 送

屿枫固位钉(膜钉)A套餐

屿枫膜钉一次购买5盒(25颗)
送
膜钉手持钉器一支(直、弯二选一)



BTA03



BTA04

或

套餐B



一次购买10盒 送



型号随机

屿枫固位钉(膜钉)B套餐
屿枫膜钉一次购买10盒(50颗)
送
滋养孔钻1.0一支
定位钻0.8一支
膜钉手持钉器(弯+直)各一支
牙骨锤一支
不可吸收膜随机型号一张



GPD08



N104



BTA03



BTA04



800/510



MAPLEISLAND MEDICAL®

HDPTFE-MEM®



Edition 2026

HDPTFE-MEM[®]



产品展示 product presentation

钛增强微孔高致密聚四氟乙烯膜
TITANIUM REINFORCED MICRO PORE HD-PTFE MEMBRANE

HDPM




HDPM1224A
12mm x 24mm



HDPM1318A
13mm x 18mm



HDPM1319A
13mm x 19mm



HDPM1424A
14mm x 24mm



HDPM1520A
15mm x 20mm



HDPM1725A
17mm x 25mm



HDPM2025A
20mm x 25mm



HDPM2030A
20mm x 30mm



HDPM2530A
25mm x 30mm



HDPM2536A
25mm x 36mm



HDPM3041A
30mm x 41mm



HDPM3040A
30mm x 40mm



HDPM3040SA
30mm x 40mm




HDPM4040SA
40mm x 40mm




HDPM4050SA
40mm x 50mm




HDPM4060SA
40mm x 60mm




HDPM5060SA
50mm x 60mm



HDPM5070SA
50mm x 70mm



HDPM6070SA
60mm x 70mm




HDPM6080SA
60mm x 80mm

HDPMB



HDPMB2021
20mm x 21mm




HDPMB2025
20mm x 25mm



HDPMB2028
20mm x 28mm



HDPMB2029
20mm x 29mm




HDPMB2127
21mm x 27mm



HDPMB2325
23mm x 25mm



HDPMB2530
25mm x 30mm




HDPMB2830
28mm x 30mm


HDPA



HDPA1525
15mm x 25mm



HDPA2035
20mm x 35mm



HDPA3045
30mm x 45mm

不可吸收微孔高致密聚四氟乙烯膜
NON-RESORBABLE MICRO PORE HD-PTFE MEMBRANE

HDPN



HDPN1520(No Core)
15mm x 20mm




HDPN2025(No Core)
20mm x 25mm



HDPN2029(No Core)
20mm x 29mm



HDPN2530(No Core)
25mm x 30mm




HDPN3040(No Core)
30mm x 40mm

HDPTFE-MEM[®]




产品特性与优势

Product Features & Advantages




高致密的聚四氟乙烯膜，
高密度的孔径防止细菌渗透

Highly dense PTFE membrane with high-density pores to prevent bacterial penetration.




本产品可以修剪，根据所需修剪形状，
但尽量不要修剪到钛支架

The membrane can be trimmed as your requirement, but try not to trim to titanium bracket.



本产品可配合钛钉、膜钉使用固定，HDP A、
HDPMB型号可配合基台螺丝使用，适配度高

The membrane can be fixed with screw or tack.
The HDP A and HDPMB models can be used with abutment screws with high compatibility.



内含钛支架可塑性强、支撑力强，可根据所需
支撑贴合

The built-in titanium scaffold has strong plasticity and strong supporting force, and can be supported and fitted according to the required.

产品注意事项

Product cautions

- 1、本产品是高致密的HDPTFE-MEM[®]
- 2、为了促进牙龈快速愈合可以在HDPTFE-MEM[®]上覆盖可吸收生物膜，但不是必选项，需结合临床自行选择
- 3、本产品不分正反面，建议有纹路的一面朝向软组织、光滑面朝向术区。产品可以用手预弯，与术区贴合

1. This is highly dense HDPTFE-MEM[®].
2. An absorbable biomembrane can be covered on the HDPTFE-MEM[®] to promote rapidly healing of gums, but it is not a must. It can be selected according to clinical situation.
3. The product has no distinction between front and back. It is recommended to place the textured side face soft tissue, and the smooth side face surgical site. It can be pre-bent by hand to conform to the surgical site.





HDPTFE-MEM[®]






产品临床案例

Product clinical cases


CASE1



CASE2



CASE3





MAPLEISLAND MEDICAL®